



# The SOLUS project: Smart optical and ultrasound diagnostics of breast cancer

**Paola Taroni**

*Politecnico di Milano, Dipartimento di Fisica*

**and the SOLUS Consortium**



## Background (I)

- **Breast cancer**

- \* The most common cancer in women in Europe and worldwide (more than 2 million new cases diagnosed in 2020) *[GLOBOCAN 2020]*
- \* About 1 in 8 women in Europe and US will be diagnosed with breast cancer in the course of her lifetime

- **Role of early diagnosis**

- \* Mammographic screening reduces breast-cancer mortality: ranging from 28% to 35% *[Nelson et al, Ann Intern Med, 2016]*

## Background (II)

- **Limitations of screening mammography**

- \* Detection of previously occult benign breast lesions
- \* The cumulative risk of a false-positive mammogram over a 10-year period of yearly screening reaches 50-60%

- **Consequences**

- \* Needless additional imaging and invasive procedures
  - Negative impact on the patient's quality of life
  - High burden for the healthcare systems

## General goal of the project

To support the *in vivo* diagnosis of breast cancer **improving the discrimination of lesions that are borderline between benign and malignant lesions**

- ➡ Patients' quality of life
- ➡ Sustainability of the healthcare systems

## “Technical” objectives towards the general goal

- **Main objective**

Develop an innovative **multi-modal tomographic system** combining *ultrasound and optical imaging*


⇒ Breast imaging    ⇒ Medical imaging

- **Specific/intermediate objective**

Develop an innovative, ground-breaking, low-cost, small size photonic module – the **Smart Optode** - to perform in depth diffuse optical measurements



⇒ Non diagnostic medical applications  
Other fields

## SOLUS concept




## Time domain diffuse optical spectroscopy

- Effects of light propagation
  - \* Attenuation
  - \* Delay
  - \* Broadening




- Theoretical model:
  - Diffusion approximation to Radiative Transport
  - Absorption coefficient → Tissue **composition**
  - Scattering coefficient → Tissue **structure**

## Absorption spectra and tissue composition






## Scattering spectra and tissue structure



## Multi-modal tomographic probe



## Diffuse Optics in SOLUS

- **Time domain**




- ❖ To uncouple absorption from scattering
- Time domain detector

- **Multiple wavelengths**

- ❖ To estimate Hb, HbO<sub>2</sub>, water, lipids, collagen
- Broad spectral range (635-1064 nm)

- **Tomographic approach**

- ❖ To identify to region of interest and improve quantification
- Time-gated wide-area detector



## State of the art of time domain diffuse optics

**2015:** 600-1350 nm  
single point




$m^3$

**2017:** 8 wavelengths  
2 channels



$dm^3$


**2020:** 8 wavelengths  
tomographic




$cm^3$

DEVICE SIZE

## Timeline



- Picosecond pulsed laser driver
- Wide area gated SiPM detector
- Dedicated acquisition electronics



**Optode components**

**Single optode**

**Multimodal probe & system**



**Clinical validation**

Developed and validated  
in laboratory settings

Integrated and  
under characterization


Start in  
Summer 2021

## Multiple wavelength picosecond pulsed lasers



- **Integrated dual-channel Laser Driver**  
Tunable pulse-width and delay with 1 ps steps  
Tested up to 80 MHz repetition frequency
- **8 wavelength lasers**  
Specifically selected (635-1064 nm)
- **Resulting performances**  
FWHM < 240 ps  
Output power: between 1.5 and 6 mW  
Negligible exponential tail or secondary peaks

## Fast-gated dSiPM with integrated TDC



- **Digital SiPM**

- Photon collection area:  $4.9 \times 4.7 \text{ mm}^2$
- Controllable active area: up to  $8.6 \text{ mm}^2$   
hot-pixel shutdown → lower noise  
signal equalization
- Gate-on transition:  $<500 \text{ ps}$ , up  $\approx 3 \text{ mm}^2$
- Temporal response  
 $235 \text{ ps}$  (FWHM) for a single active pixel  
 $\approx 500 \text{ ps}$  for an active area of  $\approx 4 \text{ mm}^2$

- **Time to Digital Converter**

78 ps resolution, 128-channel histogram builder

## The smart optode

### Ultra-compact system for time domain multi-wavelength DOT

- **Responsivity:**  $>10^6$  m<sup>2</sup>sr @600 nm for an active area of 2 mm<sup>2</sup>
  - ✓ 10x larger than for state-of-the-art systems
- Retrieval of **optical properties** (homogeneous medium)
  - ( $\mu_a = 0.06-0.4$  cm<sup>-1</sup> and reduced scattering  $\mu_s' = 4-17$  cm<sup>-1</sup>)
    - ✓ Relative error on Absorption: 10% (median)
    - Relative error on Scattering: up to 40%
- Sensitivity to **deep absorption perturbations**
  - $\Delta\mu_a = 0.16$  cm<sup>-1</sup>, 1 cm<sup>3</sup>, embedded in bkg medium ( $\mu_a = 0.1$  cm<sup>-1</sup> and  $\mu_s' = 10$  cm<sup>-1</sup>)
  - ✓ Sensitivity: down to 3.5 cm with 2% contrast



## The multi-modal SOLUS system

Calibration and performance assessment on-going



SOLUS screen


SOLUS probe



cooling system

3D position bay

medical grade power supply



Mach 30 screen

3D position sensor

standard US probe

Mach 30 body

## Next steps

### Clinical validation of the SOLUS system

- 3 radiologists
- “Mock sessions”
  - \* Train medical doctors
  - \* Test the system usability and ergonomics
- Pilot clinical feasibility study
  - \* 20 benign (BIRADS 2-3) and 20 malignant (BIRADS 4-5) lesions

### Smart optode

- Investigate the optode potential
  - \* Development of wearable devices (medical, sport, rehabilitation fields)
  - \* Use in other fields (e.g., fruit quality assessment)

# Thank you for your attention

Paola Taroni - [paola.taroni@polimi.it](mailto:paola.taroni@polimi.it)



This project has received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement No 731877. The project is an initiative of the Photonics Public Private Partnership.

[www.solus-project.eu](http://www.solus-project.eu)

[www.photonics21.org](http://www.photonics21.org)