

Correlative Light and X-ray Microscopy: CLXM

Liz Duke
Diamond Light Source
UK



Acknowledgements



Maria Harkiolaki, Kyle Dent, Matt Spink



Ben Hornberger, Dean Hunt, Jeff Irwin,
Pete Lander, Vadim Palshin, Andrei
Tkachuk



Lucy Collinson, Minoos Razi, Raffa Carzaniga,
Marie-Charlotte Domart
London Research Institute CRUK



Beamline staff on U41-TXM




Eva Pereiro


Talk Outline


- Introduction to cryo soft X-ray tomography
- Sample Preparation and Workflow
- Some results
- Going forward


A bit of background


- New 3D technique in structural biology
- Biological imaging of cells with soft X-rays
- Complementary to EM and light microscopy












Analysing Cellular Features

N: Nucleus

C: Cytoplasm



Cell membrane

Cytoskeleton

Endoplasmic reticulum

Mitochondria

Lysosomes

Nuclear envelope

Nucleoli

N

C

Where can I collect data?





Cryo-soft X-ray tomography

WHOLE CELLS

NO CHEMICAL FIXATIVES

NO STAINS

NO DEHYDRATION

'Near Native-State Imaging'

At the home laboratory


Live

- Apply for beamtime at the synchrotron
- Prepare or purchase holey-carbon grids
- Grow cells on grids and transfect if required
- Image the cells by live fluorescence microscopy

Frozen

- Plunge freeze the grids in liquid ethane and store in liquid nitrogen
- Image the frozen grids by cryo-fluorescence microscopy, noting the position of cells of interest
- Ship the grids to the synchrotron in a dry nitrogen shipper



At the synchrotron

- 
- Load grids into the microscope
 - Trace back the position of cells of interest and image using the in-vacuum visible light microscope
 - Acquire an X-ray mosaic image of the area of interest
 - Determine the final field of view and perform a test rotation to establish the extreme tilt angles and increment steps
 - Calculate dose and exposure time
 - Run the tilt series

At the home laboratory



- Import the data into e.g. IMOD software
- Align the tilt series (using fiducial markers)
- Perform BP and SIRT reconstructions
- Align the fluorescence and X-ray datasets
- Identify potential structures of interest
- Manually segment the structures and render into a 3D model
- Overlay fluorescence image with 3D model to back-check correlation

Using cryo-fluorescence imaging to screen grids prior to shipping






Van Driel, 2009

Correlative Example 1 – Fluorescence Search (20X Objective)




Correlative Example 1 – Soft X-ray Imaging




- ~4 μm thick peri-nuclear region with biological event of interest

Cryo fl using Cryostage @ LRI



In microscope VLM @ BESSY





Going Forward

- New technique:
 - Instrumentation
 - Sample Preparation
- 2D into 3D
- Super resolution